

HOF BEAUTY

Fact Sheet

1062 SkinMate Galvanic Machine




Description

This SkinMate Galvanic Machine is a great addition to any professional beauty business. This machine uses a direct galvanic mild electrical current. It is recommended for different skin types and conditions as it has scope to treat in two different ways - Desincrustation and Iontophoresis.

Desincrustation removes excess sebum debris and dead skin cells to provide a deep cleanse.

Iontophoresis allows water-soluble substances to enter unbroken skin.



HOF BEAUTY

Product Specifications

| | |
|--------------|--------------|
| Voltage | 85-265V |
| Frequency | 50-60Hz |
| Rated Power | 20W |
| Dimensions | 22*19*11.5cm |
| Package Size | 31*27*16cm |
| Net Weight | 3.4kg |

Contents

| | |
|--------------|----|
| Machine | x1 |
| Power Cable | x1 |
| Roller Probe | x1 |
| Round Probe | x1 |
| Manual | x1 |

Related Products



1075 SkinMate Digital Beauty Tower

Compatible with any SkinMate machine



3132 SkinMate Desincrustation Lotion and 3133 SkinMate Iontophoresis Lotion provide a comfortable interface during the flow of current